ENGINEERING APPLICATION 24/10/16
OF
NANOMATERIALS

AIM:

To determine the aneroge crystallite size from the given X-ray siffraction (XRD) pattern of a polycrystalline moterial.

APPARATUS REQUIRED:

XRD pattern (uphood in the course page)

Pear fitting program (Open fource / free softwar

like fityk, guplot and

Utiplot is preferrable)

is strictly prohibited.

FORMULA:

The scheme is used to calculate crystalline size here,

D = K.
Berso

DATA GIVEN:

Instrumental brooking = 0.01 M = 1.546Å sherrer noustantly 0.94

OBSERVATION TABLE:

Paol	FWMM	Fumm often instru- mental broduring correction	FUMM 2	Any oyslal sizali)
28.57		0.224 = 0.00390	0.00195	848.57
47.54	0.276	0.266=0.0046	0.00232	
56.37	0.298	0.288 = 0.00 502	0,00251	10.45.41
33.13	0.248	0.238 = 0.00415	0.002095	1836.31
59,12	0.3.5	0.295=0.00514	0.00257	(102.01

CALCULATIONS:

$$0_1 = \frac{0.94 \times 1.546}{0.00195 \times Cos(28.57)} = 848.57 \text{ Å}$$

$$D_2 = \frac{0.94 \times 1.546}{0.00232 \times \cos(47.94)} = 927.89 \text{ Å}$$

$$D_3 = \frac{0.94 \times 1.546}{0.00251 \times (200)(156.37)} = 1045.4 \text{ Å}$$

Average =
$$\frac{0.+0_2+0_3+0_4+0_5}{5}$$
 = 952.048 Å

RESULT:
The average crystalline size by scherrer formula of poly-vytalline ving
is 952.098 Å