

25/7/16

SPECTROMETER

EXP-2

AIM :

To determine the refractive index of the glass prism using a spectrometer

APPARATUS REQUIRED :

Spectrometer, spirit level, magnifying glass, glass prism, sodium vapour lamp

FORMULA :


$$\mu = \frac{\sin \left( \frac{A + \delta_m}{2} \right)}{\sin \left( \frac{A}{2} \right)}$$

where  $\mu$  = Refractive Index

$A$  = Angle of prism

$\delta_m$  = Angle of minimum deviation

DIAGRAM :





# OBSERVATION TABLE

$$\text{Least count} = \frac{\text{Length of 1 division of MS}}{\text{No. of divisions, vernier scale}} = \frac{30'}{30} = 1'$$

Angle of prism, A (as obtained earlier) =  $60^\circ$

Vernier	Reading for minimum deviation position ( $R_1$ )			Reading for direct ray ( $R_2$ )			$\delta_m = R_1 - R_2$	$\mu$
A	MSR (deg)	VSR (min)	TR (min)	MSR (deg)	VSR (min)	TR (min)	38.54'	1.51
	162°	0'	162°0'	123°	6'	123°6'		
B	MSR (deg)	VSR (min)	TR (min)	MSR (deg)	VSR (min)	TR (min)	39.2'	1.52
	342°	4'	342°4'	303°	2'	303°2'		

$$\text{Average } \mu = \frac{1.51 + 1.52}{2} = \frac{3.03}{2} = 1.515$$

CALCULATIONS:

$$\begin{aligned} A: \mu &= \frac{\sin\left(\frac{60^\circ + 38^\circ 54'}{2}\right)}{\sin(60^\circ/2)} \\ &= 2 \sin\left(\frac{98^\circ 54'}{2}\right) \\ &= 1.51 \end{aligned}$$



$$B: \mu = \frac{\sin\left(\frac{60 + 39^\circ 2'}{2}\right)}{\sin\left(\frac{60}{2}\right)}$$

$$\sin\left(\frac{60}{2}\right)$$

$$= 2 \sin(49.51)$$

$$= 1.52$$

l.c  
denver

RESULT :

The Refractive Index of the given glass prism is : 1.515

9/11/10

*[Signature]*